

MaculaPlex® 2

DHA to support normal vision



PRODUCT DESCRIPTION

Give your eyes a moment

Colours, shapes and impressions not only enrich our lives, they also provide us with important information about the world in which we live. For us humans, the eyes are the most important sensory organ. We capture 80% of our impressions with our eyes and thus create an accurate image of our outside world. However, as we age, reading, watching television or even simply recognising faces can become increasingly difficult. This is not only due to the normal aging process, but also to so-called 'free radicals', which can be created by UV rays, for example, and lead to oxidative stress. For this reason, it is essential to support and maintain our vision in the best possible way. Biogena MaculaPlex® was developed with this in mind and contains micronutrients that support the eyes and can maintain vision. It consists of two perfectly matched individual preparations and is based on the latest scientific findings: MaculaPlex® 1 contains a combination of micronutrients with an antioxidant effect, which are complemented by the omega-3 fatty acids contained in MaculaPlex® 2, and can thus make a comprehensive contribution to maintaining normal vision.

- Macula Complex with branded raw materials
- OmniXan® RR-ZEAXANTHIN & Lutemax® Free Lutein
- With active ubiquinol (KANEKA UBIQUINOL™)
- To support the eyes and vision
- With the omega-3 fatty acid DHA

POSSIBLE USES

To maintain normal eyesight, especially in old age

INGREDIENTS	PER DAILY DOSE (1 CAPSULE)	% NRV *
Fish oil concentrate	500 mg	-
contains: Docosahexaenoic acid (DHA)	250 mg	-
Vitamin E (mg alpha-TE)	8 mg	67

* % of the reference quantity according to EU Regulation 1169/2011.

INGREDIENTS

fish oil concentrate, capsule shell (gelatine from **fish**, humectant: glycerol), D-alpha-tocopherol, mixed tocopherols.

Food supplement: The recommended daily dosage should not be exceeded. Not a substitute for a balanced and varied diet and a healthy lifestyle.

Not an organic product according to VO (EU) 2018/848.

Suitable for diabetics.

STORAGE INSTRUCTIONS

Store in a dark, dry place at room temperature. Protect from heat.

Keep out of the reach of children.

QUALITY SEALS



Biogena GmbH & Co KG
Strubergasse 24, A-5020 Salzburg, biogena.com

Infoline Austria toll-free: T +43 662 23 11 11, info@biogena.com
Infoline Germany, Italy and International: T +49 8654 774 00-0, deutschland@biogena.com
Infoline Switzerland: T +41 43 215 20 59, schweiz@biogena.com

Online orders: biogena.com

